

# DICOM CONFORMANCE STATEMENT

Medical Exposure Management Information System

onti

©PDRadiopharma Inc. 2023  
ALL RIGHTS RESERVED

## Contents

|          |   |           |
|----------|---|-----------|
| <b>1</b> | <b>INTRODUCTION</b>                       | <b>2</b>  |
| 1.1      | DOCUMENT HISTORY . . . . .                | 2         |
| 1.2      | IMPORTANT . . . . .                       | 2         |
| <b>2</b> | <b>CONFORMANCE STATEMENT OVERVIEW</b>     | <b>4</b>  |
| <b>3</b> | <b>NETWORK</b>                            | <b>5</b>  |
| 3.1      | DICOM DATA FLOW DIAGRAM . . . . .         | 5         |
| 3.1.1    | Functional Definition of AE's . . . . .   | 5         |
| 3.2      | AE SPECIFICATIONS . . . . .               | 7         |
| 3.2.1    | Verification SCP . . . . .                | 7         |
| 3.2.2    | Storage SCP . . . . .                     | 9         |
| 3.2.3    | Query/Retrieve SCP . . . . .              | 12        |
| 3.2.4    | MWM SCU . . . . .                         | 15        |
| 3.3      | COMMUNICATION PROFILES . . . . .          | 17        |
| 3.4      | CONFIGURATION . . . . .                   | 17        |
| <b>4</b> | <b>SUPPORT OF EXTENDED CHARACTER SETS</b> | <b>18</b> |

# 1 INTRODUCTION

## 1.1 DOCUMENT HISTORY

| Date                   | Revision | Author             | Description     |
|------------------------|----------|--------------------|-----------------|
| October 31, 2023, 2023 | 1.0      | PDRadiopharma Inc. | Initial Version |

## 1.2 IMPORTANT

- The copyright of this document belongs to PDRadiopharma Inc.
- This content is subject to change without notice.
- Reproduction, duplication, modification, etc. of part or all of this document without permission from our company is prohibited.
- onti is a trademark of PDRadiopharma Inc. This manual may include trademarks or registered trademarks of other companies.

This document is written for the people that need to understand how the onti Software will integrate into their healthcare facility. This includes both those responsible for overall imaging network policy and architecture, as well as integrators who need to have a detailed understanding of the DICOM features of the product. This document contains some basic DICOM definitions so that any reader may understand how this product implements DICOM features. However, integrators are expected to fully understand all the DICOM terminology, how the tables in this document relate to the product's functionality, and how that functionality integrates with other devices that support compatible DICOM features. Also note that this document is formatted according to the DICOM 3.0 Specification, Part 2: Conformance.

The scope of this DICOM Conformance Statement is to facilitate integration between the "onti" Software and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard. DICOM by itself does not guarantee interoperability. The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information.

- The comparison of different Conformance Statements is just the first step towards assessing interconnectivity and interoperability between the product and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.
- The "onti" product participate in an industry-wide testing program sponsored by Integrating the Healthcare Enterprise (IHE). Together with the IHE Technical Framework, may facilitate the process of validation testing.

## 2 CONFORMANCE STATEMENT OVERVIEW

The product name is "onti". This product is a system that implements DICOM services necessary for recording and managing medical exposures. It is also possible to download the worklist from the information system. Table 1 overview the network services supported by onti.


Table 1: Network Services

| SOP Classes  | User of Service (SCU) | Provider of Service (SCP) |
|--|-----------------------|---------------------------|
| <b>Transfer</b>  |                       |                           |
| Computed Radiography Image Storage                       | No                    | Yes                       |
| Digital X-Ray Image Storage – For Presentation           | No                    | Yes                       |
| Digital Mammography X-Ray Image Storage For Presentation | No                    | Yes                       |
| CT Image Storage   | No                    | Yes                       |
| MR Image Storage   | No                    | Yes                       |
| Ultrasound Image Storage                                 | No                    | Yes                       |
| Secondary Capture Image Storage                          | No                    | Yes                       |
| X-Ray Angiographic Image Storage                         | No                    | Yes                       |
| X-Ray Radiofluoroscopic Image Storage                    | No                    | Yes                       |
| Nuclear Medicine Image Storage                           | No                    | Yes                       |
| X-Ray Radiation Dose SR Storage                          | No                    | Yes                       |
| Radiopharmaceutical Radiation Dose SR Storage            | No                    | Yes                       |
| Positron Emission Tomography Image Storage               | No                    | Yes                       |
| RT Structure Set Storage                                 | No                    | Yes                       |
| <b>Query/Retrieve</b>                                    |                       |                           |
| Study Root Information Model – FIND                      | No                    | Yes                       |
| Study Root Information Model – MOVE                      | No                    | Yes                       |
| <b>Workflow Management</b>                               |                       |                           |
| Modality Worklist  | Yes                   | No                        |

### 3 NETWORK

#### 3.1 DICOM DATA FLOW DIAGRAM

Figure 1 : DICOM DATA FLOW DIAGRAM



##### 3.1.1 Functional Definition of AE's

- Verification SCP
  - After setting the conditions of the remote AE (AE Title, IP Address, Port Number), establish an association with the remote AE using the C-ECHO protocol.
- Storage SCP
  - After setting the conditions of the remote AE (AE Title, IP Address, Port Number), an object is received from the remote AE using the C-STORE protocol.
- Query/Retrieve SCP
  - After setting the conditions of the remote AE (AE Title, IP Address, Port Number), it returns search results to the request from the remote AE using the C-FIND protocol.

- Transmits the requested image to the remote AE using the C-MOVE protocol.
- MWM SCU
  - After setting the conditions of the remote AE (AE Title, IP Address, Port Number), retrieve the Worklist from the remote AE using the C-FIND protocol.

Sequencing of Real-World Activities: Not applicable.

## 3.2 AE SPECIFICATIONS

### 3.2.1 Verification SCP

- SOP Class
  - Table2, The Verification SCP AE provides Standard Conformance to the following SOP Class.
- Association Policies
  - General
    - \* Table3, The DICOM standard application context name for DICOM 3.0 is always proposed.
  - Number of Associations
    - \* Does not support simultaneous associations. Maximum one.
  - Asynchronous Nature
    - \* Does not support asynchronous communication.
- Association Acceptance Policy
  - Negotiates for an Association establishment request by C-ECHO of the remote AE (SCU).
  - Accepts a request from a remote AE (SCU) to establish an association by C-ECHO and returns a "normal" or "error" response to the remote AE (SCU) after the internal processing.
- Accepted Presentation Contexts
  - Table4, Accept the Presentation Context.
- SOP specific conformance for Verification SOP Class
  - Provides standard conformance to the Verification SCP Class that can be received as an SCP.
  - Use C-STORE



Table 2: Standard Conformance to the following SOP Class

| <b>SOP Class Name</b> | <b>SOP Class UID</b> | <b>SCP</b> |
|-----------------------|----------------------|------------|
| Verification          | 1.2.840.10008.1.1    | Yes        |

Table 3: Application Context Name

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

Table 4: Presentation Context

| <b>Abstract Syntax</b> |                   | <b>Transfer Syntax</b>    |                   | <b>Role</b> | <b>Extended Negotiation</b> |
|------------------------|-------------------|---------------------------|-------------------|-------------|-----------------------------|
| <b>Name</b>            | <b>UID</b>        | <b>Name</b>               | <b>UID</b>        |             |                             |
| Verification C-ECHO    | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP         | None                        |

### 3.2.2 Storage SCP

- SOP Class
  - Table5, The Storage SCP AE provides Standard Conformance to the following SOP Class.
- Association Policies
  - General
    - \* Table6, The DICOM standard Application Context Name for DICOM is always accepted.
  - Number of Associations
    - \* Does not support simultaneous associations. Maximum one.
  - Asynchronous Nature
    - \* Does not support asynchronous communication.
- Association Acceptance Policy
  - Receives a request by C-STORE of a remote AE (SCU).
  - Receives a request by C-STORE from a remote AE (SCU) for a Receivable class Object (Images and Dose SR).
  - Received objects are stored in the specified local location.
- Accepted Presentation Contexts
  - Accepts the Presentation Context of Table7 and Table8.
- SOP specific conformance for Storage SOP Class
  - Provides standard conformance to the Storage SOP Class that can be received as an SCP.
  - Use C-STORE

Table 5: Standard Conformance to the following SOP Class

| SOP Classes  | SOP Class UID                 | SCP |
|--|-------------------------------|-----|
| Computed Radiography Image Storage                         | 1.2.840.10008.5.1.4.1.1.1     | Yes |
| Digital X-Ray Image Storage – For Presentation             | 1.2.840.10008.5.1.4.1.1.1.1   | Yes |
| Digital Mammography X-Ray Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.2   | Yes |
| CT Image Storage   | 1.2.840.10008.5.1.4.1.1.2     | Yes |
| MR Image Storage   | 1.2.840.10008.5.1.4.1.1.4     | Yes |
| Ultrasound Image Storage                                   | 1.2.840.10008.5.1.4.1.1.6.1   | Yes |
| Secondary Capture Image Storage                            | 1.2.840.10008.5.1.4.1.1.7     | Yes |
| X-Ray Angiographic Image Storage                           | 1.2.840.10008.5.1.4.1.1.12.1  | Yes |
| X-Ray Radiofluoroscopic Image Storage                      | 1.2.840.10008.5.1.4.1.1.12.2  | Yes |
| Nuclear Medicine Image Storage                             | 1.2.840.10008.5.1.4.1.1.20    | Yes |
| X-Ray Radiation Dose SR Storage                            | 1.2.840.10008.5.1.4.1.1.88.67 | Yes |
| Radiopharmaceutical Radiation Dose SR Storage              | 1.2.840.10008.5.1.4.1.1.88.68 | Yes |
| Positron Emission Tomography Image Storage                 | 1.2.840.10008.5.1.4.1.1.128   | Yes |
| RT Structure Set Storage                                   | 1.2.840.10008.5.1.4.1.1.481.3 | Yes |

Table 6: Application Context Name

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

Table 7: Abstract Syntax

| Name List  | UID List                      | Role | Extended Negotiation |
|--|-------------------------------|------|----------------------|
| Computed Radiography Image Storage                         | 1.2.840.10008.5.1.4.1.1.1     | SCP  | None                 |
| Digital X-Ray Image Storage – For Presentation             | 1.2.840.10008.5.1.4.1.1.1.1   | SCP  | None                 |
| Digital Mammography X-Ray Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.2   | SCP  | None                 |
| CT Image Storage   | 1.2.840.10008.5.1.4.1.1.2     | SCP  | None                 |
| MR Image Storage   | 1.2.840.10008.5.1.4.1.1.4     | SCP  | None                 |
| Ultrasound Image Storage                                   | 1.2.840.10008.5.1.4.1.1.6.1   | SCP  | None                 |
| Secondary Capture Image Storage                            | 1.2.840.10008.5.1.4.1.1.7     | SCP  | None                 |
| X-Ray Angiographic Image Storage                           | 1.2.840.10008.5.1.4.1.1.12.1  | SCP  | None                 |
| X-Ray Radiofluoroscopic Image Storage                      | 1.2.840.10008.5.1.4.1.1.12.2  | SCP  | None                 |
| Nuclear Medicine Image Storage                             | 1.2.840.10008.5.1.4.1.1.20    | SCP  | None                 |
| X-Ray Radiation Dose SR Storage                            | 1.2.840.10008.5.1.4.1.1.88.67 | SCP  | None                 |
| Radiopharmaceutical Radiation Dose SR Storage              | 1.2.840.10008.5.1.4.1.1.88.68 | SCP  | None                 |
| Positron Emission Tomography Image Storage                 | 1.2.840.10008.5.1.4.1.1.128   | SCP  | None                 |
| RT Structure Set Storage                                   | 1.2.840.10008.5.1.4.1.1.481.3 | SCP  | None                 |

Table 8: Transfer Syntax

| <b>Name List</b>  | <b>UID List</b>        |
|---|------------------------|
| Implicit VR – Little Endian                             | 1.2.840.10008.1.2      |
| Explicit VR – Little Endian                             | 1.2.840.10008.1.2.1    |
| JPEG Lossless, Nonhierarchical, First- Order Prediction | 1.2.840.10008.1.2.4.70 |
| JPEG-LS Lossless Image Compression                      | 1.2.840.10008.1.2.4.80 |
| JPEG 2000 Image Compression (Lossless)                  | 1.2.840.10008.1.2.4.90 |
| JPEG 2000 Image Compression (Lossy)                     | 1.2.840.10008.1.2.4.91 |

### 3.2.3 Query/Retrieve SCP

- SOP Class
  - Table9, The Query/Retrieve SCP AE provides Standard Conformance to the following SOP Class.
- Association Policies
  - General
    - \* Table10, The DICOM standard Application Context Name for DICOM is always accepted.
  - Number of Associations
    - \* Does not support simultaneous associations. Maximum one.
  - Asynchronous Nature
    - \* Does not support asynchronous communication.
- Association Acceptance Policy
  - Returns matching information for a C-FIND request from a remote AE (SCU).
  - When a C-MOVE request is received from a remote AE (SCU), the result of the request is returned.
  - For a C-FIND request from a remote AE (SCU), the contents of the request are checked against the internal database and the matching information is returned.
  - When a C-MOVE request is received from a remote AE (SCU), it returns the result (success/failure) of the transfer of the request contents.
- Accepted Presentation Contexts
  - Accepts the Presentation Context of Table11.
- SOP specific conformance for Query SOP Class
  - 3 levels (Study, Series, Image) are supported for C-FIND requests. The standard search item is Table12.
- SOP specific conformance for Retrieve SOP Class
  - Use C-MOVE
  - The transfer syntax supported Table13.

Table 9: Standard Conformance to the following SOP Class

| SOP Classes                         | SOP Class UID               | SCP |
|-------------------------------------|-----------------------------|-----|
| Study Root Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Yes |
| Study Root Information Model – MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Yes |

Table 10: Application Context Name

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

Table 11: Presentation Context

| Abstract Syntax                     |                             | Transfer Syntax           |                   | Role | Extended Negotiation |
|-------------------------------------|-----------------------------|---------------------------|-------------------|------|----------------------|
| Name                                | UID List                    | Name                      | UID List          |      |                      |
| Verification                        | 1.2.840.10008.1.1           | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP  | None                 |
| Study Root Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP  | None                 |
| Study Root Information Model – MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP  | None                 |

Table 12: Search

| Level  | Attribute Name     | Tag        | Type         |
|--------|--------------------|------------|--------------|
| Study  | Study Date         | 0008, 0020 | Required Key |
|        | Study Time         | 0008, 0030 | Required Key |
|        | Accession Number   | 0008, 0050 | Required Key |
|        | Study ID           | 0020, 0010 | Required Key |
|        | Study Instance UID | 0020, 000D | Unique Key   |
|        | Study Description  | 0008, 1030 | Option       |
| Series | Modality           | 0008, 0060 | Required Key |
|        | Series Number      | 0008, 0011 | Required Key |
| Image  | SOP Instance UID   | 0008, 0018 | Unique Key   |
|        | Instance Number    | 0020, 0013 | Required Key |

Table 13: Transfer Syntax

| <b>Name List</b>   | <b>UID List</b>        |
|--|------------------------|
| Implicit VR – Little Endian                              | 1.2.840.10008.1.2      |
| Explicit VR – Little Endian                              | 1.2.840.10008.1.2.1    |
| JPEG Lossless, Nonhierarchical, First – Order Prediction | 1.2.840.10008.1.2.4.70 |
| JPEG-LS Lossless Image Compression                       | 1.2.840.10008.1.2.4.80 |
| JPEG 2000 Image Compression (Lossless)                   | 1.2.840.10008.1.2.4.90 |
| JPEG 2000 Image Compression (Lossy)                      | 1.2.840.10008.1.2.4.91 |

### 3.2.4 MWM SCU

- SOP Class
  - Table14, The MWM SCU AE provides Standard Conformance to the following SOP Class.
- Association Policies
  - General
    - \* Table15, The DICOM standard Application Context Name for DICOM is always accepted.
  - Number of Association
    - \* Does not support simultaneous associations.
    - \* One Association is started at a time for each worklist request.
  - Asynchronous Nature
    - \* Does not support asynchronous communication.
- Association Acceptance Policy
  - A Worklist query request is initiated.
  - Requests C-FIND. Send to remote AE (SCP).
  - Get response from the remote AE (SCP).
- Accepted Presentation Contexts
  - Accepts the Presentation Context of Table16.
- SOP specific conformance for Modality Worklist
  - Attributes are DICOM standard.



Table 14: Standard Conformance to the following SOP Class

| <b>SOP Class</b>                           | <b>SOP Class UID</b>   | <b>SCU</b> |
|--|------------------------|------------|
| Modality Worklist Information Model – FIND | 1.2.840.10008.5.1.4.31 | Yes        |

Table 15: Application Context Name

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

Table 16: Presentation Context

| <b>Abstract Syntax</b>                     |                        | <b>Transfer Syntax</b>    |                   | <b>Role</b> | <b>Extended Negotiation</b> |
|--|------------------------|---------------------------|-------------------|-------------|-----------------------------|
| <b>Name</b>                                | <b>UID</b>             | <b>Name</b>               | <b>UID</b>        |             |                             |
| Modality Worklist Information Model – FIND | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU         | None                        |

### **3.3 COMMUNICATION PROFILES**

- Supported Communication Stack
  - DICOM Part 8: TCP/IP Network Communication
- TCP/IP Stack
  - Use the TCP/IP stack of the Microsoft Windows Operating System.

### **3.4 CONFIGURATION**

- Mapping
  - AE Title, TCP/IP and Port Number are set.
- Parameter
  - Setting
    - \* AE Title (Default: onti\_ae)
    - \* TCP/IP
    - \* Port Number
    - \* MWM – SCU (Search acquisition item settings.)

## 4 SUPPORT OF EXTENDED CHARACTER SETS

Table17, This product supports the following character sets.

Table 17: character sets

| <b>ISO</b>      |                    |
|-----------------|--------------------|
| ISO 2022 IR 6   | Default (ISO 646)  |
| ISO 2022 IR 87  | JIS X 0208 (Kanji) |
| ISO 2022 IR 100 |                    |