

Techne® MDP Kit

Indications / Effects	<ul style="list-style-type: none">•Diagnosis of bone diseases by bone scintigraphy•Diagnosis of cerebral tumor or cerebral vessel disorders by cerebral scintigraphy
Package	2 vials, 10 vials (lyophilized)
Contents	One vial contains: MDP (methylene diphosphonate): 7.5 mg Tin (II) chloride dehydrate: 0.38 mg Ascorbic acid: 0.17 mg Hydrochloric acid: Appropriate amount Sodium hydroxide: Appropriate amount
pH	5.0 – 6.0 (after dissolution by Sodium pertechnetate (^{99m} Tc) injection)
Storage	2 – 8°C
Expiration	24 months after production