

IHE Integration Statements

IHE Integration Statement			
Vendor	Product Name	Version	Date
PDRadiopharma Inc.	Bridgea GATEWAY	1.0.0 or later	2022-12-27
This product implements all transactions required in the IHE Technical Framework to support the IHE profiles, actors and options listed below.			
Profiles Implemented	Actors Implemented	Connectathon Seal	Options Implemented
Radiation Exposure Monitoring for Nuclear Medicine (REM-NM)	Radiopharmaceutical Activity Supplier (RAS)		None
Consistent Time (CT)	Time Client		None
Internet address for vendor's IHE information: https://www.pdradiopharma.com/products/5009/			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	https://www.pdradiopharma.com/products/5009/		
Links to general information about IHE worldwide: http://ihe.net/IHE Worldwide			