







# IHE Integration Statements

IHE Integration Statement			
Vendor	Product Name	Version	Date
PDRadiopharma Inc.	ankan	1.0.0 or later	2022-12-21
This product implements all transactions required in the IHE Technical Framework to support the IHE profiles, actors and options listed below.			
Profiles Implemented	Actors Implemented	Connectathon Seal	Options Implemented
Radiation Exposure Monitoring (REM)	Dose Information Consumer		None
	Image Manager / Archive		None
Radiation Exposure Monitoring for Nuclear Medicine (REM-NM)	Dose Information Consumer		None
	Dose Information Reporter		None
	Image Manager / Archive		None
Consistent Time (CT)	Time Client		None
<b>Internet address for vendor's IHE information:</b> <a href="https://www.pdradiopharma.com/products/6456/">https://www.pdradiopharma.com/products/6456/</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	None		
<b>DICOM</b>	<a href="https://www.pdradiopharma.com/products/6456/">https://www.pdradiopharma.com/products/6456/</a>		
<b>Links to general information about IHE worldwide:</b> <a href="http://ihe.net/IHE_Worldwide">http://ihe.net/IHE_Worldwide</a>			